



68.02.08
17,5 cm – 7"
Obwegeser



68.02.09
7 mm

68.02.10
11 mm




68.02.00
17,5 cm – 7"
Hourigan



für die Ramus-Vorderkante
for the ascending ramus
para el borde delantero del ramus
pour le bord antérieur du rameau
per lo spigolo del ramo


68.02.05
19 cm – 7½"
Hargis






68.26.15
17 cm – 6¾"
Steinhäuser

für den Orbitaboden
for the orbital floor
para la base de la orbita
pour le fond de l'orbite
per il pavimento orbitale



68.26.12

18 cm – 7⅛"
Tessier



68.26.13



für den linken Unterkieferrand
for left inferior border of the mandible
para el borde izquierdo de la mandíbula
pour bord gauche du maxillaire inférieur
per il bordo sinistro della mandibola

68.02.06
18 cm – 7⅛"
Obwegeser



für den rechten Unterkieferrand
for right inferior border of the mandible
para el borde derecho de la mandíbula
pour bord droit du maxillaire inférieur
per il bordo destro della mandibola

68.02.07
18 cm – 7⅛"
Obwegeser



68.02.21
 21 cm – 8¼"
Obwegeser-Freer




68.02.22
 20,5 cm – 8"
Obwegeser-Freer




15 cm – 6"
Chompret




329.80.01
 # 1



329.80.02
 # 2




329.80.03
 # 3



329.80.04
 # 4





68.03.02

2,5 mm

68.03.04

4,0 mm

68.03.39

6,5 mm

15 cm – 6"

Obwegeser

Interdentale Osteotome
Interdental osteotomes
Osteótomos interdentaes
Ostéotomes interdentaux
Osteotomi interdentali



7 mm

68.03.07

21,5 cm – 8½"

Obwegeser

Spaltungsosteotome
Splitting osteotomes
Osteótomos hendedores
Ostéotomes à fendre
Osteotomi fenditori



10 mm

68.03.10

21,5 cm – 8½"

Obwegeser



10 mm

68.03.30

23,5 cm – 9¼"

Obwegeser

Tuber Osteotom
Pterygoid osteotome
Osteótomo pterigoideo
Ostéotome pterygoïde
Osteotomo per tubercolo



**Osteotome
Osteotomes
Osteótomos
Ostéotomes
Osteotomi**




68.03.12
8 mm

68.03.13
12 mm

68.03.14
16 mm

22 cm – 8¾"
Obwegeser

Keilosteotome
Wedge osteotomes
Osteótomos de cuña
Ostéotomes cunéiformes
Osteotomi a forma di cuneo



68.03.16
8 mm

68.03.17
12 mm

68.03.18
16 mm

68.03.19
20 mm

22 cm – 8¾"
Obwegeser

Blattosteotome
Blade osteotomes
Osteótomos de hoja
Ostéotomes en lame
Osteotomi a forma di lama



68.03.34
8 mm

68.03.35
10 mm

21 cm – 8¼"
Obwegeser

flexibel
flexible
flexible
flexible
flessibile




68.26.17
5 mm

68.26.18
8 mm

22,5 cm – 9"
Steinhäuser

Orbital Osteotome
Orbital osteotomes
Osteótomos orbitales
Ostéotomes orbitaux
Osteotomi orbitali






68.03.50
10 mm

68.03.52
15 mm

68.03.54
20 mm

68.03.56
30 mm

20,5 cm – 8"
Tessier



68.03.51
10 mm

68.03.53
15 mm

20,5 cm – 8"
Tessier



68.03.60
10 mm

68.03.62
15 mm

16 cm – 6¼"
Tessier



68.03.61
10 mm


68.03.63
15 mm

16 cm – 6¼"
Tessier



68.03.75

Nasenseptum Osteotome
Nasal septum osteotomes
Osteótomos para el tabique nasal
Ostéotomes pour cloison nasale
Osteotomi per setto nasale



68.03.76

18,5 cm – 7¼"
Obwegeser



Repositionszangen, Fragment-Adaptierungsklemmen
Reposition Forceps, Fragment Adapting Forceps
Pinzas para reposicionar, Pinzas adaptadoras para fragmentos
Pinces à réposition, Clamps d'adaptation de fragments osseux
Pinze per riposizione, Pinze per avvicinare frammenti



35.11.13
 13,5 cm – 5 $\frac{3}{8}$ "





35.11.14
 14,5 cm – 5 $\frac{3}{4}$ "



68.25.88
 13 cm – 5 $\frac{1}{8}$ "
Steinhäuser




68.25.89
 17,5 cm – 7"
Steinhäuser




68.40.02
 19,5 cm – 7 $\frac{3}{4}$ "
Obwegeser


rechts
 right
 derecha
 droite
 destra





68.40.01
 19,5 cm – 7 $\frac{3}{4}$ "
Obwegeser

links
 left
 izquierda
 gauche
 sinistra






68.40.04



68.40.05

22 cm – 8¾"
Obwegeser



68.40.10
23,5 cm – 9¼"
Smith




Positionierer für Knochenfragmente, Mobilisator
Positioner for Bone Fragments, Mobilizer
Posicionador de fragmentos, Mobilizador
Positionneur des fragments osseux, Mobilisateur
Posizionatore di frammenti ossei, Mobilizzatore

Positionierer für UK-Knochenfragmente
 Positioner for mandibular bone fragments
 Posicionador de fragmentos mandibulares
 Positionneur des fragments osseux mandibulaires
 Posizionatore di frammenti ossei della mandibola



68.05.20
 18 cm – 7 1/8"
Wangerin

Maxilla-Mobilisator
 Maxillary mobilizer
 Mobilizador maxilar
 Mobilisateur pour maxillaire
 Mobilizzatore per mascella



68.03.29
 23 cm – 9"
Obwegeser





68.04.05
 14 cm – 5½"
Obwegeser

für Unterkiefer
 for mandibles
 para mandíbulas
 pour maxillaires inférieures
 per mandibole



68.04.10
 23 cm – 9"
Obwegeser

für Jochbögen
 for zygomatic archs
 para arcos cigomáticos
 pour arcades cygomatique
 per zigomi



68.04.15
 18 cm – 7⅛"
Obwegeser

für Oberkiefer
 for maxillas
 para maxilares
 pour maxillaires supérieures
 per mascelle




Dingman

10 mm



77.10.10
 18 cm – 7⅛"


20 mm



77.10.20
 17 cm – 6¾"




**Elevatorien
Elevators
Elevadores
Elévateurs
Elevatori**



scharf
sharp
cortante
tranchant
tagliente

stumpf
blunt
roma
mousse
smusso


26.37.12
12 S
17 cm – 6¾"
Bennet



scharf
sharp
cortante
tranchant
tagliente

scharf
sharp
cortante
tranchante
tagliente

26.37.42
22
19 cm – 7½"
Seldin




scharf
sharp
cortante
tranchant
tagliente


stumpf
blunt
roma
mousse
smusso

26.37.43
23
19,5 cm – 7¾"
Seldin






68.09.05
 18 cm – 7 1/8"
Dingman




68.45.10
 23 cm – 9"
Gillies



68.09.20
 21 cm – 8 1/4"
Dingman




für Unterkiefer
for mandibles
para mandíbulas
pour maxillaires inférieures
per mandibole







Repositionszangen
Reposition Forceps
Pinzas para reposicionar
Pinces à réposition
Pinze per riposizione




68.50.01
 24 cm – 9½"




68.50.02
 24 cm – 9½"



68.50.03
 23,5 cm – 9¼"



68.50.04
 23,5 cm – 9¼"



68.50.05
 22,5 cm – 9"




Rowe

links
 left
 izquierda
 gauche
 sinistra






Rowe

rechts
 right
 derecha
 droite
 destra




68.50.06
 22,5 cm – 9"






Raspatorien
 Raspatories
 Legras
 Rugines
 Periostotomi







Lippen- und Gaumenspalten-Instrumente
Lip- and Cleft Palate Instruments
Instrumentos para labio y fisuras paladares
Instruments pour lèvre et fissures palatines
Strumenti per labbra e fessura del palato



68.82.79
 15,5 cm – 6 1/8"
Wang




68.82.80
 15,5 cm – 6 1/8"
Wang



rechts
 right
 derecha
 droite
 destra



68.82.81
 18 cm – 7 1/8"
Wang



links
 left
 izquierda
 gauche
 sinistra



68.82.82
 18 cm – 7 1/8"
Wang



rechts
 right
 derecha
 droite
 destra



68.82.83
 18 cm – 7 1/8"
Wang



links
 left
 izquierda
 gauche
 sinistra



68.82.84
 18 cm – 7 1/8"
Wang

Raspatorien
 Raspatories
 Legras
 Rugines
 Periostotomi





68.82.35
 # 1
 19 cm – 7½"



68.82.36
 # 2
 19 cm – 7½"




Muko-Periostal Haken
 Mucoperiost retractors
 Separadores mucoperiostales
 Ecarteurs mucoperiosteals
 Divaricatori mucoperiostali


68.82.37
 # 3
 19 cm – 7½"



Fräser, Tamponstopfer, Drahtinstrumente
Burrs, Gauze Packer, Wire Instruments
Fresas, Tamponador, Instrumentos para alambre
Fraises, Porte-tampons, Instruments pour fils métalliques
Frese, Portatampone, Strumenti per fili d'acciaio




Lindemann



Drahtdrehklemme
Wire twister
Pinza para torcer alambre
Pince à torsader le fil
Pinza per torcere fili

Drahtösenklemme
Wire loop forceps
Pinza para ojete de alambre
Pince pour œillets en fil métallique
Pinza per occhiello in filo metallico





**68.15.05
Erich**

Rolle/1 m
Roll/1 m
Rollo/1 m
Rouleau/1 m
Rot/1 m

Drahtbögen
Arch bars
Arcos de alambre
Arcs en fil
Archi di filo



68.15.10
14 cm – 5½"
Erich

Verpackungseinheit 12 Stück
Packing unit 12 pieces
Unidad de 12 piezas
Emballage de 12 pièces
Confezione da 12 pezzi



Drahtbügelschiene
Arch bar
Arco de alambre
Rail en fil
Arco di filo


68.04.93
15 cm – 6"
Dautrey



76.01.05
Ø 0,5 mm

Knochendraht, weich
Bone wire, soft
Alambre de huesos, dulce
Fil à os, souple
Filo per ossa, dolce

Rolle/10 m
Roll/10 m
Rollo/10 m
Rouleau/10 m
Rot/10 m



**68.15.01
Schuchardt**

Verpackungseinheit 10 Stück
Packing unit 10 pieces
Unidad de 10 piezas
Emballage de 10 pièces
Confezione da 10 pezzi

Sprossenschienen
Lattice splints
Rieles en forma de rejilla
Rails à échelons
Rotaie a forma di pinolo




**68.15.02
Schuchardt**


Verpackungseinheit 10 Stück
Packing unit 10 pieces
Unidad de 10 piezas
Emballage de 10 pièces
Confezione da 10 pezzi



Beckenkamm-Retraktoren
Iliac Crest Retractors
Separadores para cresta iliaca
Ecarteurs iliaques
Divaricatori per cresta iliaca



68.01.67
30/60 × 60 mm






68.01.68
30/60 × 80 mm



22 cm – 8¾"
Obwegeser








Sinuslifting Elevatorien
Sinuslifting Elevators
Elevadores para sinus-lift
Elévateurs pour sinus-lift
Elevatori per sinus-lifting


Raspatorium
 Raspatory
 Legar
 Rugine
 Periostotomo




26.37.08
 18 cm – 7 1/8"
Stiller




20 cm – 8"
Stiller






68.82.53
 3 mm




68.82.54
 4 mm



68.82.55
 6 mm



Komprimierer
 Compressor
 Compresor
 Compresseur
 Compessore





Zylinder-Osteotome
Cylinder Osteotomes
Osteotomos cilíndricos
Ostéotomes cylindriques
Osteotomi cilindrici

Kortikalisstanze
 Cortical cutter
 Punzonadora cortical
 Perforateur cortical
 Perforatore corticale



Ø 12 mm



78.41.02
 12 cm - 4¾"
Shepard


Querstift
 Cross pin
 Espiga transversal
 Embout transversal
 Perno trasversale



zum Drehen der Zylinder-Osteotome
 for turning cylinder osteotomes
 para girar osteótomos cilíndricos
 pour la rotation de l'ostéotome cylindrique
 per la rotazione degli osteotomi cilindrici



78.41.11
 11 cm - 4⅜"
Shepard



Zylinder-Osteotome
 Cylinder osteotomes
 Osteotomos cilíndricos
 Ostéotomes cylindriques
 Osteotomi cilindrici



Shepard



12,5 cm - 5"

Ø = 5 mm **78.41.05**

Ø = 6 mm **78.41.06**

Ø = 8 mm **78.41.08**

Ausstoßer
 Ram
 Expulsor
 Pousoir
 Espulsore



14 cm - 5½"

78.41.12

78.41.13

78.41.14

Trokare
 Trocars
 Trocares
 Trocart
 Trocar



14,5 cm - 5¾"

78.41.15

78.41.16

78.41.18



Knochenmühle mit Montageschlüssel
Bone mill with assembling wrench
Trituradora de hueso con llave de montaje
Moulin à os avec clé de montage
Trita ossa con chiave per montaggio



Ermöglicht ein Zerkleinern von Knochenstücken
Enables crushing of compact bone pieces
Permite moler fragmentos óseos
Pour broyage des fragments d'os
Rende possibile tritare frammenti ossei