



2.4


Titan-Schrauben für MANDIBULAR-System mit Werkzeug
Titanium Screws for MANDIBULAR System with Tools
Tornillos de titanio para sistema MANDIBULAR con sus respectivos instrumentos
Vis en titane pour système MANDIBULAR avec outils
Viti in titanio per sistema MANDIBULAR con strumentario

“POWER GRIP”

Schraubenkopf mit Kreuzschlitz
Screw head with cross slot
Tornillo de cabeza ranurada en cruz
Tête de vis avec fente cruciforme
Testa della vite con taglio a croce




P₁₀



Verpackungseinheit 10 Schrauben im Magazin (nicht steril).
Packing unit 10 screws in storage tray (non-sterile).
Paquete de 10 tornillos en cartucho (no estéril).
Emballage unitaire de 10 vis en cartouche (non stériles).
Confezione da 10 viti nella cartuccia (non sterile).

Spiralbohrer
Twist drills
Brocas espirales
Forets hélicoïdaux
Perforatori elicoidali



Ansatz
Shank end
Talón
Embout
Attacco



STRYKER


7 - 23	7 - 23	7 - 23
a = 25 mm b = 50 mm 68.95.72	a = 25 mm b = 80 mm 68.95.64	a = 25 mm b = 105 mm 68.95.69



Verpackungseinheit “X” Stück
Packing unit “X” pieces
Paquete de “X” unidades
Emballage de “X” pièces
Confezione da “X” pezzi



Selbstschneidend
Self-tapping
Autoroscante
Autotaraudeuse
Autofilettante




	Self-tapping  Ø 2,4 mm	Self-tapping with locking  Ø 2,4 mm	Emergency screw  Ø 2,7 mm
c			
7,0 mm	68.94.50A	68.94.80A	68.94.60A
9,0 mm	68.94.51A	68.94.81A	68.94.61A
11,0 mm	68.94.52A	68.94.82A	68.94.62A
13,0 mm	68.94.53A	68.94.83A	68.94.63A
15,0 mm	68.94.54A	68.94.84A	68.94.64A
17,0 mm	68.94.55A	68.94.85A	68.94.65A
19,0 mm	68.94.56A	68.94.86A	
21,0 mm	68.94.57A	68.94.87A	
23,0 mm	68.94.58A	68.94.88A	

Schraubendreher
Screwdrivers
Destornilladores
Tournevis
Cacciavite



Selbstschneidend mit Verblockung
Self-tapping with locking
Autoroscante con blocco
Autotaraudeuse avec blocage
Autofilettante con bloccaggio



Notfallschraube (selbstschneidend)
Emergency screw (self-tapping)
Tornillo de emergencia (autoroscante)
Vis d'urgence (autotaraudeuse)
Vite d'emergenza (autofilettante)